

PROPRIETARY MACHINE VISION IDENTIFICATION ALGORITHM



AnCam™ 6100 MULTI-ANALYTE TEST READER

Features

- Reads most single and multi-target coupons
- Machine vision minimizes human error
- Archives results, time and GPS location
- Transmits data and images via local cellular network
- 24-hour continuous assays
- Gives early warning of strong positives
- Detects defective coupons, incorrectly inserted coupons, or no control line
- Water resistant, durable case
- Compact design

Application Areas

- Field detection of bioterrorism agents
- Medical pathogens
- Drugs of abuse
- Food safety
- Veterinary testing
- Environmental testing

THE ANCAM™ 6100 is an extremely versatile portable multi-analyte test reader and cellular communication device usable with virtually any single or multi-target lateral flow coupon. Lateral flow immunoassay tests called coupons are used to detect and identify bioterrorism threat agents, medical pathogens, drugs of abuse, environmental toxins, food-borne pathogens, and animal disease, to name a few. Results are determined by visual inspection of the coupon. However, human vision is imperfect, particularly in situations of high stress or low light. A simple visual examination may also not detect manufacturing defects, out-of-date reagents, and errors in procedure. The AnCam 6100 replaces the human eye with an independent evaluation of the entire test protocol as well as the results. Human and ambient lighting issues are removed.

The AnCam 6100 Multi-Analyte Test Reader uses a modified smart phone to provide a highly integrated solution that combines machine vision, proprietary signal processing algorithms, GPS-based locating capabilities, data storage and data transmission in one compact package, along with standard smart phone capabilities such as texting, email, and Bluetooth-enabled communications. This smart phone-based system can also detect signals not easily visible to the eye, increasing test sensitivity. With calibration, quantitative measurements can be obtained instead of a binary Yes/No result.

The AnCam 6100 helps the user perform each test correctly and provides preliminary estimates as the test progresses. For example,

many tests of this type are assumed to require a 15-minute incubation period. But the AnCam's proprietary machine vision algorithms provide unprecedented sensitivity and reliability. In most cases where there is a reasonable concentration of the target material, a positive reading will become apparent as early as five minutes into the incubation, reducing the threat detection response time by a factor of three.

The AnCam 6100 can be used to analyze coupons that have been previously incubated, or for monitoring the development process in a new coupon from the point of fluid sample introduction to final outcome. Fluid sample flows and the concentration of colored reagents within all coupon channels are monitored to determine that the proper amount of fluid sample was used and that the coupon incubation protocol is proceeding normally.

Once the assay incubation period is complete, results may be transmitted to a remote location as a text message or email using the local cellular network. All images and assay calculations involving the assay are stored in the AnCam 6100's internal memory for later examination, and optionally on its SD card.

AnCam 6100 Specifications

| | |
|------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------|
| Lateral flow coupon compatibility | Most commercially available tests including 1, 5, and 8 analyte coupons |
| Image recognition/processing | Proprietary high resolution algorithms |
| Types of analysis | Single measurement or tracked incubation Qualitative or quantitative results |
| Protocol identification | Recognition of visual coupon features |
| Camera and machine vision engine | Customized embedded cellular phone provides GPS/map functions |
| Processing time | Results within seconds |
| Archival capabilities | Date and time, user ID, coupon information, GPS latitude/longitude, one or more coupon photographs saved for later analysis |
| Communication | Wi-Fi and cellular network; send assay result message to selected cell phone numbers or email addresses |
| Hardcopy output | Via Bluetooth printer or Wi-Fi printer |
| Size | 205 x 121 x 109mm |
| Weight | 1280 grams |
| Temperature range | 0° to 50°C |
| Power | Internal lithium-ion battery |
| Standby operating time | 5-10 days |

CBRN International, Ltd. reserves the right to change specifications without prior notice.

For more information or to request a price quote, please contact us at the email address below.

CBRN International
Rigga Business Center, Unit 5001
Al Rigga Road, Deira
Dubai P.O. Box 4647
Dubai, UAE

info@cbrnintl.com
1.425.922-4007
Registration No.:
186350 Jafza Offshore

